## KMS LOW AIR LOSS /ALTERNATING PRESSURE MATTRESS SYSTEM

WITH DELUXE PUMP AND CELL-ON-CELL SUPPORT BASE




- Dual therapies; low air loss and alternating pressure
- Unique cell-on-cell design ensures that mattress will hold air in lower cells in the event of power failure
- Nylon top cover is waterproof and quilted for added comfort
- Weight capacity 450 lb



1-800-488-6535

## KMS Low Air Loss Alternating Pressure Mattress



## Bladders slowly deflate and inflate to constantly change and re-distribute pressure points

Alternating Pressure Mattresses can help heal or prevent decubitus ulcers, commonly called "bedsores" or "pressure ulcers" and provide comfort to patients with burns, Multiple Sclerosis, Lou Gehrig's disease, Fibromyalgia, impaired mobility etc.

KMS low air loss mattresses have cell-on-cell design and bladders with tiny pin-holes which diffuse air up into the space between the cells and the vapor-permeable top cover. Low air loss mattresses float the patient on air-filled cells while circulating air across the skin to reduce moisture and to help maintain a constant, skin interface pressure. Low air loss can reduce the negative effects of heat and moisture buildup on the skin.

KMS King Medical Supply

1-800-488-6535